

Protein (Urine)

Interpretive Summary

Description: Protein in the urine can be an indicator of renal disease (interpret in conjunction with urine concentration) or lower urinary tract disease, when pre-renal causes of proteinuria have been ruled out.

Negative Protein

Common Causes

- Clinically normal animal
- False negatives:
 - Bence Jones protein unreliable
 - Highly buffered alkaline urine

Increased Protein

Common Causes

- Small amounts of protein (50 mg/dL or less) can be normal in urine, especially if urine is concentrated
 - Interpret results in conjunction with urine specific gravity
 - Further evaluation by urine protein:creatinine ratio to determine clinical significance of observed proteinuria.
- Prerenal proteinuria
 - Overflow/overload pre-glomerular proteinuria
 - Bence Jones proteins, hyperglobulinemia
 - Post-colostral proteinuria
 - Hemoglobinuria, myoglobinuria
 - Systemic hypertension
- Renal proteinuria
 - Glomerular proteinuria
 - Glomerulonephritis
 - Amyloidosis
 - Tubular proteinuria
 - Acute renal disease
 - Defects in proximal renal tubular function
 - Congenital disorders
 - Functional – transient
 - Fever, shock, cardiac disease, exercise (horse), seizures
- Postrenal Proteinuria
 - Inflammation or infection of the upper or lower urinary tract (or reproductive tract in voided specimen)
 - Hemorrhage into the urinary tract or reproductive tract
 - Neoplasia of the urinary tract
- Cushing's disease
- False positives
 - Dipstick
 - Prolonged contact of dipstick reagent pad with alkaline urine
 - Quaternary ammonium or chlorhexidine contamination
 - Pigmenturia
 - Myoglobin
 - Hemoglobin
 - Recent administration of certain blood substitutes (e.g. Oxyglobin)
 - Sulfosalicylic acid test
 - Radiographic contrast media
 - High doses of antibiotics (penicillin, cephaloridin, sulfisoxazole)

- Co-precipitation of urinary crystals due to low pH of reagent

Related Findings

- Prerenal proteinuria
 - Increased globulins
 - Positive Bence Jones proteins in urine
 - Hemoglobinuria
 - Anemia
 - Increased bilirubin
 - Myoglobinuria
 - Increased CK, AST
- Renal proteinuria
 - Increased urine protein:creatinine ratio
 - Decreased albumin
 - Secondary systemic hypertension
 - Increased BUN, creatinine, phosphorus with secondary renal tubular damage
 - Increased cholesterol, ascites/pulmonary edema in severe nephrotic syndrome
 - Positive serologic/PCR testing if glomerulonephritis secondary to infectious agents
 - Heartworm, Lyme, leptospirosis, rickettsial, fungal, protozoal, or viral (in some cases)
 - Positive antinuclear antibody titer if associated with systemic immune-mediated disease
 - Consistent renal biopsy and electron microscopy results
- Postrenal proteinuria
 - Active urine sediment exam (RBC, WBC, bacteria, abnormal epithelial cells)
 - Positive urine culture
 - Positive bladder tumor analyte test (canine only)
 - Urinary calculi visualized on radiographs or ultrasound
- Cushing's Disease
 - Increased ALP
 - Decreased urine specific gravity
 - Adrenal function tests consistent with Cushing's disease

Additional Information

Diagnostic Methodology

- Semiquantitative
 - Colorimetric/dipstick (reagent strip)
 - Used as a screening test
 - Primarily detects albumin; does not reliably detect globulins or Bence Jones proteins associated with multiple myeloma.
 - Reported as negative, trace, 1+ to 3+ reaction that correlates to 100, 300, or 500 mg/dL protein
 - Acid precipitation tests (includes sulfosalicylic acid test - SSA)
 - Detects albumin and nonalbumin proteins, including Bence Jones paraprotein
 - Commonly used to confirm dipstick results
- Quantitative
 - Colorimetric, spectrophotometric
 - Electrophoresis, immunoelectrophoresis

References

- Feldman EC, Nelson RW. *Canine and Feline Endocrinology and Reproduction*, 3rd ed. St. Louis, MO: Saunders; 2004.
- Kaneko JJ, Harvey JW, Bruss ML. *Clinical Biochemistry of Domestic Animals*, 6th ed. San Diego, CA: Academic Press; 2008.

- Latimer KS, Mahaffey EA, Prasse KW, eds. *Duncan and Prasse's Veterinary Laboratory Medicine: Clinical Pathology*, 4th ed. Ames, IA: Blackwell; 2003.
- Osborne CA, Stevens JB. *Urinalysis: A Clinical Guide to Compassionate Patient Care*. Shawnee Mission, KS: Bayer Corporation; 1999.
- Stockham SL, Scott MA. *Fundamentals of Veterinary Clinical Pathology*, 2nd ed. Ames, IA: Blackwell; 2008.
- Willard MD, Tvedten H, eds. *Small Animal Clinical Diagnosis by Laboratory Methods*, 4th ed. St. Louis, MO: Saunders; 2004.

Last updated 11/1/2013