Catalyst One* Chemistry Analyzer

Reference Guide





Navigating the Analyzer



Analyzer Status

Description
Analyzer is ready to process samples or perform maintenance tasks
Analyzer is in standby mode
Analyzer is processing a sample or performing another activity
Analyzer is waiting for the user to begin processing a sample (after receiving the patient information from the IDEXX VetLab* Station)
An error has occurred (review alert messages on the IDEXX VetLab Station)

Note: You can also view the analyzer status by viewing its icon on the IDEXX VetLab Station Home screen.

Compatible Species

Canine [†]	Lizard	Rabbit	Avian (budgerigar, cockatoos
Feline [†]	Llama	Rat	[grey cheek, Moluccan, umbrella], cockatiel, canary,
Equine [†]	Monkey	Sea Turtles	conure, macaw
Bovine (beef cattle, dairy cow)	Mouse	Sheep	[blue and gold, hyacinth,
Ferret	Other	Snake	scarlet], parrots [Amazon blue, Amazon yellow, eclectus,
Goat	Pig	Tortoise	African grey])

[†]Species- and life stage-specific intervals are available for these species.

Analyzing a Sample

Prepare the sample. For detailed instructions, see the Preparing Samples For Use on IDEXX Catalyst* Chemistry Analyzers poster.

Sample Volume Sample Type Required Whole blood $0.7 \text{ cc} (700 \,\mu\text{L})$ Fill to the lowest line on the separator Serum. Varies depending on the number of slides plasma, or urine being run. (300 μ L is the recommended sample volume when running a Chem 17 CLIP and electrolytes.)

Initiate the run on the IDEXX VetLab Station.



Load the sample, slides, pipette tips, and any other necessary materials in the sample drawer, as shown above.

Note: Do not load a whole blood separator AND a sample cup for a single run.



Close the sample drawer and press the Start button on the analyzer.




Advanced Testing and Other Work Flows

For detailed instructions on running quality control, cleaning the internal components of the analyzer monthly, or running automated dilutions and UPC, see the on-screen instructions provided on the IDEXX VetLab Station touch screen.

Catalyst CLIPs and Slides

Chemistry	Abbreviation	Chem 17 CLIP	Chem 15 CLIP	Chem 10 CLIP	Equine 15 CLIP	NSAID 6 CLIP	UPC Panel⁺	Lyte 4 CLIP	SDMA and TT ₄ Kit	QC CLIP	Individual Slides
Albumin	ALB	✓	✓	✓	✓					✓	✓
Alkaline phosphatase	ALKP	✓	✓	✓	✓	✓				✓	✓
Alanine aminotransferase	ALT	✓	✓	✓		✓				✓	✓
Amylase	AMYL	✓									✓
Aspartate aminotransferase	AST				✓	✓					✓
Bile acids	BA										✓
Blood urea nitrogen	BUN	✓	✓	✓	✓	✓					✓
Calcium	Ca	✓	✓		✓					✓	✓
Cholesterol	CHOL	✓	✓								✓
Creatine kinase	СК				✓						✓
Creatinine	CREA	✓	✓	✓	✓	✓					✓
Chloride	CI							✓			
C-reactive protein [‡]	CRP										✓
Fructosamine [†]	FRU										✓
Gamma-glutamyltransferase	GGT	✓	✓		✓						✓
Glucose	GLU	✓	✓	✓	✓					✓	✓
Potassium	К							✓			
Lactate	LAC										✓
Lactate dehydrogenase	LDH				✓						✓
Lipase	LIPA	✓									✓
Magnesium	Mg										✓
Sodium	Na							✓			
Ammonia	NH ₃									✓	✓
Phenobarbital [†]	PHBR										✓
Inorganic phosphate	PHOS	✓	✓								✓
Progesterone	PROG										✓
Symmetric dimethylarginine [†]	SDMA								✓		✓
Total bilirubin	TBIL	✓	✓		✓						✓
Total protein	TP	✓	✓	✓	✓						✓
Total T ₄ [†]	TT4								✓		✓
Triglycerides	TRIG										✓
Urine creatinine	UCRE						✓				
Urine protein	UPRO						✓				
Uric acid	URIC										✓



[‡]Validated reference intervals for feline, equine, and "other" species are unavailable.

