**SNAP® Combo Plus**  
 English version

SNAP Combo Plus FeLV Ag/FIV Ab test is an enzyme-linked immunosorbent assay (ELISA) for the simultaneous detection of feline leukemia virus (FeLV) antigen and antibodies to feline immunodeficiency virus (FIV) in feline serum, plasma or whole blood.

- Precautions and warnings.**
- All waste should be properly decontaminated prior to disposal.
  - Do not mix components from kits with different serial numbers.
  - Do not use a SNAP device that has been activated prior to the addition of sample.
  - Refer to country specific Material Safety Data Sheet for regional hazard identification.

**Storage**

- SNAP devices and test reagents must be stored at 2-8°C.
- All components must be at room temperature 18-25°C before running the test - **DO NOT HEAT.**

**Kit components**

Item	Reagent	5 test Quantity	15 test Quantity	30 test Quantity
1	1 bottle Anti-FelV/FIV Ag HRPO conjugate	3.5 mL	3.5 mL	7.0 mL
2	SNAP device: Monoclonal antibodies to p27, inactivated FIV antigen, and positive and negative controls	5	15	30

**Reagents contained in each device:**

Wash solution	0.4 mL	0.4 mL	0.4 mL
Substrate solution	0.6 mL	0.6 mL	0.6 mL

**Other components:** transfer pipettes, sample tubes, and reagent rack

**WARNING**

Conjugate - H317/H412/P280/P302 + P332/P333 + P313: May cause an allergic skin reaction. Harmful to aquatic life with long lasting effects. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

Wash Solution - H317/P280/P302 + P332/P333 + P313: May cause an allergic skin reaction. Wear protective gloves. IF ON SKIN: Wash with plenty of soap and water. If skin irritation or rash occurs: Get medical advice/attention.

- Sample information.**
- Samples must be at room temperature 18-25°C before beginning test procedure.
  - Serum or plasma, either fresh, previously frozen or stored at 2-8°C, may be used in this test.
  - Serum or plasma may be stored for up to 7 days at 2-8°C. For longer storage, sample should be frozen (-20°C or colder).
  - Previously frozen or older samples must be centrifuged before use.
  - Hemolyzed or lipemic samples will not affect results.
  - EDTA or heparin in plasma will not affect results.
  - Whole blood may be used. Whole blood must be anticoagulated (eg. EDTA, heparin) and may be used either fresh or after refrigeration 2-8°C for up to one week.

**Test procedure**

- Using the pipette that is provided, transfer **3 drops of sample** (whole blood, serum, or plasma) into sample tube.
  - Holding bottle vertical, add **4 drops of conjugate** to the sample tube.
  - Cap the sample tube and mix thoroughly by **inverting tube 3-5 times**.
  - Place the device on a flat surface. Add contents of sample tube to Sample Well, being careful not to splash contents outside of Sample Well. Sample will flow across Result Window, reaching Activator in approximately 30-60 seconds. Some sample may remain in Sample Well. **Watch carefully for sample or blue color in the Activate Circle.**
  - When color **FIRST** appears in Activate Circle push Activator firmly until it is flush with the device body.
- 
- Note:** Some samples may not flow to the activate circle within 60 seconds, and, therefore, the circle may not turn color. In this case, press the Activator after sample has flowed across Result Window.
- Keep the device horizontal to ensure accurate results.
  - Read test result at **10 minutes**. Note: positive control may develop sooner, but results are not complete until 10 minutes.

**Interpreting test result**

To determine test results, read the reaction spots in the Result Window. Color development in sample spots is proportional to the concentration of FeLV antigen or FIV antibody in the sample. If no color develops in the positive control spot, repeat the test.

**Negative result**  
 Only positive control spot develops color.

**Positive result**  
 FeLV Antigen  
 Positive control spot and FeLV Ag sample spot develop color.

**FeLV antigen and FIV antibody**  
 Positive control spot and both sample spots develop color.

**FIV antibody**  
 Positive control spot and FIV Ab sample spot develop color.

**Reaction with negative control**  
 The negative control spot serves as a safeguard against false positives and helps indicate that the assay has been run properly.


**Positive result**  
 If color in the FIV Ab or FeLV Ag sample spot is darker than negative control spot, result is positive for that spot.

- Invalid results**
- Background**  
 If the sample is allowed to flow past the Activate Circle, background color may result. Some background color is normal. However, if colored background obscures test result, repeat the test.
  - No color development**  
 If positive control does not develop color, repeat the test.

**IDEXX SNAPShot Dx® Analyzer**  
 Test results can also be read using the SNAPShot Dx Analyzer. A complete description of how to enter patient data and read test results using the SNAPShot Dx Analyzer can be found in the SNAPShot Dx Analyzer Operator's Guide.

**Symbol descriptions**

Symbol	Use by date	EC REP	Authorized representative in the European Community
[LOT]	Batch code (lot)	[E]	
[SN]	Serial number	[I]	Consult instructions for use
[T]	Temperature limitation	[VD]	In vitro diagnostic
[M]	Manufacturer	[D]	Date of manufacture
[REF]	Catalog number		

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