IDEXX VetTrol* Control

For in vitro use

For veterinary use only.

Product Description

VetTrol Control is a lyophilized, bovine-based control fluid, specially formulated for use on IDEXX VetTest*, VetLyte* and Catalyst* analyzers.

Intended Use

VetTrol Control is designed for use in monitoring the precision and accuracy of the VetTest, VetLyte and Catalyst analyzers. After reconstitution, the control serum should be tested in the same manner as a patient sample, and the reported values compared to those on the printout.

Transportation

VetTrol Control is not classified as a dangerous substance for transportation purposes. However, provide adequate protection to prevent breakage and ship under refrigerated conditions.

Storage

Control and diluent vials should be stored frozen (-18°C). Discard vials at their expiration dates. Expired or unwanted material should be disposed of with other clinical waste.

Stability and Handling

For most constituents, VetTrol Control can be used up to 24 hours after reconstitution when it is stored at 2°–8°C (36°–46°F). For creatine kinase and ammonia values, VetTrol should be used within two hours after reconstitution. Exposure to light will affect bilirubin and creatine kinase results. Ammonia concentration will increase with time.

Quality Control Procedure

- Remove one diluent and one control vial from the freezer, and allow 60–90 minutes for the vials to equilibrate to room temperature.
 - Slowly invert the diluent vial several times to thoroughly mix the contents. DO NOT SHAKE.
 - Gently tap the control vial on the counter several times to dislodge any material adhering to the stopper.
 - Remove the seal and stopper from each vial just before adding the diluent. Do not leave vials open.
 - Transfer EXACTLY 3.0 mL of diluent to the control vial, using a clean, dry, Class A volumetric pipette, or an equivalent automatic pipette. A 3.0 mL syringe, with the needle removed, can be used, but the reduced volumetric precision may cause inaccurate results. Discard the remaining diluent.
 - Replace the original stopper on the control vial and hold it firmly in place. Gently invert the vial. DO NOT SHAKE.
- Reconstitution, with occasional inversion, takes 45 to 60 minutes. Visually verify that all freeze-dried material is dissolved before using.**
- Keep all control fluids tightly stoppered and refrigerated when not in use. Ensure that the fluid is at room temperature before use. Vials should be inverted several times just before using.
 - Discard the control vial after one day.

Troubleshooting

Catalyst analyzer

Refer to the Catalyst* Operator’s Manual for troubleshooting information.

VetTest analyzer/VetLyte analyzer

If your control results do not fall within the specified range for the lot of quality-control material:

- Check the expiration dates of the control material and the VetTest chemistry slides.
- Confirm that the lot number on the control vial corresponds to the lot number selected on the VetTest analyzer.
- Rerun the chemistry in question using the same control.
- If results do not fall within the specified range, rerun the chemistry using a new control vial.
- If results remain out of range, call IDEXX Technical Support or contact your authorized IDEXX representative.

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*VetTrol, VetTest, VetLyte and Catalyst are trademarks or registered trademarks of IDEXX Laboratories, Inc. in the United States and/or other countries.
All rights reserved. • 06-02248-04
© 2008 IDEXX Laboratories, Inc. Todos los derechos reservados.
IDEXX Europe B.V.
© 2008 IDEXX Laboratories, Inc. Tutti i diritti sono riservati.